



# Mouse anti-Human CGREF1 monoclonal antibody, clone 3E8 (CABT-B9957)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	CGREF1 (AAH34764, 136 a.a. ~ 235 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	3E8
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB,IP,sELISA,ELISA
<b>Sequence Similarities</b>	PAELINFPGVALRHVEPGEPLAPSPQEPQAVGRQSLAKSPLRQETQEAPGPREEAKGQV EARRESLDPVQEPGGQAEADGDVPGPRGEAEGQAEAKGDA
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In 1x PBS, pH 7.2
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Keywords** CGREF1; cell growth regulator with EF-hand domain 1; CGR11; cell growth regulator with EF hand domain protein 1; hydrophobestin; cell growth regulatory gene 11 protein; cell growth

regulator with EF hand domain 1;

---

## GENE INFORMATION

Entrez Gene ID [10669](#)

---

UniProt ID [Q99674](#)

---

Function calcium ion binding

---