



Anti-SKIV2L2 polyclonal antibody (CABT-BL033)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein of human SKIV2L2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Molecular Weight	118kDa
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Superkiller viralicidic activity 2-like 2 is a protein that in humans is encoded by the SKIV2L2 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID [23517](#)

UniProt ID [P42285](#)
