



Anti-CES2A (aa 254-303) polyclonal antibody (CABT-BL1052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: GVALLPDLIS DTSEVVYKTV ANLSGCEATD SEALIHCLRA KSKQEILAIN, corresponding to a region within internal sequence amino acids 254-303 of Mouse CES6
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	102022
Protein Refseq	NP_001177259

UniProt ID

[Q8QZR3](#)
