



Anti-CYP2D1 polyclonal antibody (CABT-BL1187)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IHC, WB
Format	Serum
Preservative	None
Storage	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction	debrisoquine 4-hydroxylase.
--------------	-----------------------------

GENE INFORMATION

Entrez Gene ID	266684
Protein Refseq	NP_695225

UniProt ID

[P10633](#)
