



Anti-DES monoclonal antibody, clone 6B7 (CABT-BL1246)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Acetic acid-solubilized desmin.
Isotype	IgA
Source/Host	Mouse
Species Reactivity	Manduca sexta
Clone	6B7
Conjugate	Unconjugated
Applications	IHC
Preservative	None
Storage	-20°C, Avoid freeze / thaw cycles

BACKGROUND

Introduction	Desmin are class-III intermediate filaments found in muscle cells. In adult striated muscle they form a fibrous network connecting myofibrils to each other and to the plasma membrane from the periphery of the Z-line structures.
--------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

GENE INFORMATION

Entrez Gene ID	395906
----------------	------------------------