



Anti-FAM59A (C-terminal) polyclonal antibody (CABT-BL1476)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to the C terminal of Human FAM59A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Buffer	99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	This gene encodes an adaptor protein which functions in the epidermal growth factor (EGF) receptor-mediated signaling pathway. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]
--------------	---

GENE INFORMATION

Entrez Gene ID [64762](#)

Protein Refseq [NP_001229338](#)

UniProt ID [Q9H706](#)
