



Anti-Gl24 (N-terminal) polyclonal antibody (CABT-BL1621)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to N-terminal residues of Gl24(Mouse)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Membrane.
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	74048
Protein Refseq	NP_001153044

UniProt ID

[E9PUF5](#)
