



Anti-GNG2 monoclonal antibody, clone EML1501 (CABT-BL1640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length protein: MASNTASIA QARKLVEQLK MEANIDRIKV SKAAADLMAY CEAHAKEDPL LTPVPASENP FREKKFFCAI L, corresponding to amino acids 1-72 of Human GNG2
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	EML1501
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-2 is a protein that in
---------------------	--

humans is encoded by the GNG2 gene.

GENE INFORMATION

Entrez Gene ID	54331
Protein Refseq	NP_001230702
UniProt ID	P59768
