



Anti-ADCK4 monoclonal antibody, clone 2E0 (CABT-BL210)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ADCK4 (AAH13114, 445 a.a. ~ 545 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2E0
Conjugate	Unconjugated
Applications	WB,ELISA
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Preservative	None
Storage	Store at -2°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene encodes a protein with two copies of a domain found in protein kinases. The encoded protein has a complete protein kinase catalytic domain, and a truncated domain that contains only the active and binding sites of the protein kinase domain, however, it is not known whether the protein has any kinase activity. Multiple transcript variants encoding different
--------------	--

isoforms have been found for this gene.

GENE INFORMATION

Entrez Gene ID	79934
----------------	-----------------------

Protein Refseq	NP_001136027
----------------	------------------------------

UniProt ID	A0A024R0Q9
------------	----------------------------
