



Anti-ALOX8 polyclonal antibody (CABT-BL258)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Mouse tissues Alox8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB, IHC, IP, ELISA
Format	Liquid
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	11688
Protein Refseq	NP_033791.1
UniProt ID	O35936