



Anti-PCNP monoclonal antibody, clone 3808D3b (CABT-BL2863)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment derived from internal sequence of human PCNP
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3808D3b
Purification	Protein G purified
Conjugate	Unconjugated
Applications	Dot, WB
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction May be involved in cell cycle regulation.

GENE INFORMATION

Entrez Gene ID [57092](#)

Protein Refseq [NP_065090](#)

UniProt ID [Q8WW12](#)
