



Anti-PLEKHF1 (aa 243-257) polyclonal antibody (CABT-BL2923)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region between amino acids 243-257 of Human APPD.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus. Cytoplasm, perinuclear region. Lysosome. Note=Translocates to lysosome during apoptosis.
Format	Liquid
Buffer	PBS, Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [79156](#)

Protein Refseq [NP_077286](#)

UniProt ID [Q96S99](#)
