



Anti-PRDM13 (aa 649-698) polyclonal antibody (CABT-BL2999)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: VKSRHPGQSL LAKAGDGPGA EPGYPPEPGD PKSDSDVDV CFTDDQSDPE, corresponding to C terminal amino acids 649-698 of Human PRDM13
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [59336](#)

Protein Refseq [NP_067633](#)

UniProt ID [Q9H4Q3](#)
