



Anti-Prolactin monoclonal antibody, clone 283 [HRP] (CABT-BL3014)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified Recombinant Human Prolactin (>98%)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	283
Conjugate	HRP
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Stable for at least one year when stored at 2°C - 8°C. Do not freeze.

BACKGROUND

Introduction	This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Aug 2011]
--------------	---

GENE INFORMATION

Entrez Gene ID [5617](#)

Protein Refseq [NP_000939](#)

UniProt ID [P01236](#)
