



Anti-SALL2 monoclonal antibody, clone 4G8 (CABT-BL3225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SALL2 (AAH24245, 1 a.a. ~ 199 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	4G8
Conjugate	Unconjugated
Applications	WB,IHC,sELISA,ELISA
Buffer	In 1x PBS, pH 7.2
Preservative	None
Storage	Store at -2°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Sal-like protein 2 is a protein that in humans is encoded by the SALL2 gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID	6297
-----------------------	----------------------

Protein Refseq [NP_005398](#)

UniProt ID [B4DK65](#)
