



Anti-SLC37A1 polyclonal antibody (CABT-BL3343)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified Human SLC37A1 recombinant protein fragment (17 kDa)
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	IgY fraction
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; multi-pass membrane protein
Format	Liquid
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	SLC37A1, a member of the sugar-phosphate transport family, transports glycerol-3-phosphate (G3P) between cellular compartments for its utilization in several compartment-specific biochemical pathways.
---------------------	--

GENE INFORMATION

Entrez Gene ID [54020](#)

Protein Refseq [NP_061837](#)

UniProt ID [P57057](#)
