



Anti-SRRM4 (aa 270-319) polyclonal antibody (CABT-BL3440)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 270-319 (KTASPLTTSR GRSQEYDSGN DTSSPPSTQT SSARSRGQEK GSPSGGLSKS) of Human SRRM4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	SRRM4 promotes alternative splicing and inclusion of neural-specific exons in target mRNAs
---------------------	--

GENE INFORMATION

Entrez Gene ID	84530
----------------	-----------------------

Protein Refseq	NP_919262
----------------	---------------------------

UniProt ID	A7MD48
------------	------------------------
