



# Anti-TAS1R2.1 (internal region) polyclonal antibody (CABT-BL3510)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide from the internal part of Zebrafish TAS1R2a conjugated to an immunogenic carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Cellular Localization</b>	Cell membrane; Multi pass membrane protein
<b>Format</b>	Liquid
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [664691](#)

Protein Refseq [NP\\_001034920](#)

---

UniProt ID [Q2MHL1](#)

---