



## Anti-TBC1D3H (aa 13-42) polyclonal antibody (CABT-BL3530)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide between 13-42 amino acids from the N-terminal region of Human TBC1D3H
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

### GENE INFORMATION

Entrez Gene ID	<a href="#">729877</a>
Protein Refseq	<a href="#">NP_001116864</a>

UniProt ID

[P0C7X1](#)

---