



Anti-TIAM1 polyclonal antibody (CABT-BL3581)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptides from human T-lymphoma invasion and metastasis-inducing protein 1 (TIAM-1, TIAM1) conjugated to an immunogenic carrier protein was used as the antigen.
Source/Host	Goat
Species Reactivity	Human
Purification	Whole serum
Conjugate	Unconjugated
Applications	IHC, WB
Format	Lyophilised
Preservative	None
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Ship	This item will be shipped to you at ambient temperature in a lyophilised form.

BACKGROUND

Introduction	Modulates the activity of RHO-like proteins and connects extracellular signals to cytoskeletal activities. Acts as a GDP-dissociation stimulator protein that stimulates the GDP-GTP exchange activity of RHO-like GTPases and activates them. Activates RAC1, CDC42, and to a lesser extent RHOA.
---------------------	--

GENE INFORMATION

Entrez Gene ID [7074](#)

Protein Refseq [NP_003244.2](#)

UniProt ID [Q13009](#)
