



Anti-ZNF236 (internal region) polyclonal antibody (CABT-BL3883)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to a region within the internal sequence 1079-1128 (VSATGETEGGDICMEEEEHSDRNASRKSRPEVITFTEE ETAQLAKIRPQ) of Human ZNF236, NP_031371
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	7776
-----------------------	----------------------

Protein Refseq [NP_031371](#)

UniProt ID [Q9UL36](#)
