



Anti-CXORF21 (full length) polyclonal antibody (CABT-BL4458)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Human CXORF21.
Immunogen	Full length human CXORF21 protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID	80231
Protein Refseq	NP_079435
UniProt ID	Q9HAI6
Chromosome Location	Xp21.3