



Anti-FGD6 (full length) polyclonal antibody (CABT-BL4663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Human FGD6.
Immunogen	Full length FGD6 protein (Human).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Purified, liquid.
Concentration	lot specific
Size	50 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID	55785
Protein Refseq	NP_060821
UniProt ID	Q6ZV73
Chromosome Location	12q23.1
Function	Rho guanyl-nucleotide exchange factor activity; guanyl-nucleotide exchange factor activity; metal ion binding; small GTPase binding;