



## Anti-FGD6 (full length) polyclonal antibody (CABT-BL4663)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Polyclonal antibody to Human FGD6.
<b>Immunogen</b>	Full length FGD6 protein (Human).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Purification</b>	Affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Purified, liquid.
<b>Concentration</b>	lot specific
<b>Size</b>	50 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

### GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">55785</a>
<b>Protein Refseq</b>	<a href="#">NP_060821</a>
<b>UniProt ID</b>	<a href="#">Q6ZV73</a>
<b>Chromosome Location</b>	12q23.1
<b>Function</b>	Rho guanyl-nucleotide exchange factor activity; guanyl-nucleotide exchange factor activity; metal ion binding; small GTPase binding;