



Anti-C11H16ORF80 (aa 144-193) polyclonal antibody (CABT-BL500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from a region between residues 144 - 193 of human GTL3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	415627
Protein Refseq	NP_001005833
UniProt ID	Q5ZHP3