



# Anti-PAMR1 (internal region) polyclonal antibody (CABT-BL5355)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit polyclonal antibody to Human PAMR1.
<b>Immunogen</b>	Synthetic peptide derived from the internal domain of Human PAMR1 protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Cellular Localization</b>	Secreted
<b>Format</b>	Liquid
<b>Size</b>	200 µl
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze/thaw cycles.

## GENE INFORMATION

Entrez Gene ID [25891](#)

<b>Protein Refseq</b>	<a href="#">NP_056245</a>
<b>UniProt ID</b>	<a href="#">Q6UXH9</a>
<b>Chromosome Location</b>	11p13
<b>Function</b>	catalytic activity; serine-type endopeptidase activity;