



# Anti-PCBD2 (full length) polyclonal antibody (CABT-BL5359)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse polyclonal antibody to Human PCBD2.
<b>Immunogen</b>	PCBD2 full-length human protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Sequence Similarities</b>	Belongs to the pterin-4-alpha-carbinolamine dehydratase family.
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.

## GENE INFORMATION

Entrez Gene ID [84105](#)

<b>Protein Refseq</b>	<a href="#">NP_115527</a>
<b>UniProt ID</b>	<a href="#">Q9H0N5</a>
<b>Chromosome Location</b>	5q31.1
<b>Function</b>	4-alpha-hydroxytetrahydrobiopterin dehydratase activity; lyase activity; phenylalanine 4-monooxygenase activity; protein binding;