



Anti-TAS2R60 polyclonal antibody (CABT-BL5818)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to Human TAS2R60.
Immunogen	A synthetic peptide as a part of human T2R60 conjugated to immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Size	100 µl
Buffer	PBS, Whole serum
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID [338398](#)

Protein Refseq	NP_803186
UniProt ID	P59551
Chromosome Location	7
Pathway	Taste transduction, organism-specific biosystem; Taste transduction, conserved biosystem;
Function	G-protein coupled receptor activity; receptor activity; signal transducer activity;