



Anti-THRSP (N-terminal) polyclonal antibody (CABT-BL5845)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to Mouse THRSP.
Specificity	Mainly expressed in tissues that synthesize triglycerides.
Immunogen	Synthetic peptide corresponding to N terminal residues of mouse THRSP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Sequence Similarities	Belongs to the SPOT14 family.
Cellular Localization	Nucleus. Cytoplasm.
Format	Liquid
Size	100 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID	21835
Protein Refseq	NP_033407
UniProt ID	Q543J4
Function	protein binding; protein homodimerization activity;