



Anti-THYMOSIN β 9 (aa 1-4) polyclonal antibody (CABT-BL5877)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	polyclonal Antibody to THYMOSIN β 9 (aa 33-41)
Specificity	Bovine Thymosin β 9 (aa 33-41)
Immunogen	Synthetic bovine Thymosin β 9 (aa 33-41) KLH conjugated (TIEQEKQAK).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Conjugate	Unconjugated
Format	Serum
Concentration	20 μ l/100 μ l (lyophilized) resuspend in 20 μ l/100 μ l aqua bidest
Preservative	None
Storage	2°C-8°C (lyophilized); -20°C (dissolved) Repeated thawing and freezing must be avoided

GENE INFORMATION

Entrez Gene ID	282385
Protein Refseq	NP_777048
UniProt ID	P21752

Functionactin binding
