



Anti-TUBB6 (aa 146-174) polyclonal antibody (CABT-BL5927)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to Human TUBB6.
Immunogen	A KLH conjugated synthetic peptide between 146-174 amino acids from the central region of Human TUBB6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Sequence Similarities	Belongs to the tubulin family.
Cellular Localization	Cytoplasm; cytoskeleton.
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	84617
Protein Refseq	NP_115914
UniProt ID	Q9BUF5
Chromosome Location	18p11.21
Pathway	Chaperonin-mediated protein folding, organism-specific biosystem; Cooperation of Prefoldin and TriC/CCTin actin and tubulin folding, organism-specific biosystem; Formation of tubulin folding intermediates by CCT/TriC, organism-specific biosystem; Gap junction, organism-specific biosystem; Gap junction, conserved biosystem; Metabolism of proteins, organism-specific biosystem; Pathogenic Escherichia coli infection, organism-specific biosystem; Pathogenic Escherichia coli infection, conserved biosy
Function	GTP binding; GTPase activity; nucleotide binding; structural molecule activity;