



Anti-ZMAT5 (aa 1-170) polyclonal antibody (CABT-BL6019)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to Human ZMAT5.
Immunogen	Recombinant Human ZMAT5 protein 1-170aa (BC009717).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Sequence Similarities	Contains 1 C3H1-type zinc finger.
Cellular Localization	Nucleus.
Format	Liquid
Size	100 µg
Buffer	pH: 7.20 99% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	55954
Protein Refseq	NP_061976
UniProt ID	A0A024R1I1
Chromosome Location	22cen-q12.3
Function	metal ion binding; nucleic acid binding; zinc ion binding;