



Anti-OYDV polyclonal antibody (CABT-BL6154)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Isotype	IgG
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live of 2 years after date of purchase.

BACKGROUND

Introduction	A virus with flexuous filamentous particles c. 775 nm long. Readily transmitted by several aphid species in the non-persistent manner and by inoculation of sap, mainly to certain Allium species. It is probably distributed world-wide.
Keywords	Onion yellow dwarf virus; OYDV; Allium virus 1; Marmor cepae