



Anti-PVA polyclonal antibody (CABT-BL6163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Isotype	IgG
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live of 2 years after date of purchase.

BACKGROUND

Introduction	An RNA-containing virus with filamentous particles, normal length c. 730 nm. Only hosts known are in the Solanaceae. Transmitted by aphids in the non-persistent manner and by inoculation of sap. Widely distributed in potato-growing areas.
Keywords	Potato virus A; PVA; Marmor solani; Potato mild mosaic virus; Potato virus P; Solanum virus 3