

This product is for research use only and is not intended for diagnostic use.

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Isotype	IgG
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live off-live ofفوق ١٠ سنوات

<b>Introduction</b>	A virus with RNA-containing tubular rod-shaped particles c. 308 x 18 nm, easily transmitted by mechanical inoculation and by handling during cultivation. No natural vector known. Widely
---------------------	---

distributed throughout the world. Of economic importance for tobacco crops.

---

**Keywords**

Tobacco mild green mosaic virus; TMGMV; Mild dark-green tobacco mosaic virus; Para-tobacco mosaic virus; Mild strain of tobacco mosaic virus; South Carolina mild mottling strain of tobacco mosaic virus; Strains U2 and U5 of tobacco mosaic virus

---