



Anti-Dd polyclonal antibody (CABT-BL6190)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using pure cultures of bacteria as immunogen.
Isotype	IgG
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA, IF
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live of 2 years after date of purchase.

BACKGROUND

Introduction	Ech (or Dickeya dianthicola, the species that has been found in GB) has the potential to cause high levels of wilting in potato crops. This is particularly evident in the warmer southerly potato growing regions of Europe. There is also evidence that low tuber inoculum levels can cause significant disease. The factors influencing disease development are generally the same as for blackleg.
---------------------	--

Keywords

Dickeya dianthicola (Erwinia chrysanthemi); Dd
