



Anti-Xhc polyclonal antibody (CABT-BL6211)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using pure cultures of bacteria as immunogen.
Isotype	IgG
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.1 M phosphate buffered saline, pH 7.4.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live off-live of 6 months.

BACKGROUND

Introduction	Bacterial blight of carrot (<i>Daucus carota</i>) is caused by <i>Xanthomonas hortorum</i> pv. <i>carotae</i> . The pathogen is seed transmitted and carrot seeds can be an important source of primary inoculum.
Keywords	<i>Xanthomonas hortorum</i> pv <i>carotae</i> ; Xhc; <i>Xanthomonas campestris</i> pv. <i>Carotae</i> ; <i>Xanthomonas carotae</i>