



# Anti-CymMV polyclonal antibody [AP] (CABT-BL6237)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using purified virus as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

## BACKGROUND

Introduction	Cymbidium mosaic virus and the Odontoglossum ringspot virus (ORSV) are two of the most common viruses affecting cultivated orchids worldwide. Infected plants can have less desirable flowers or other problems, causing significant financial losses to orchid growers. The virus has not often been reported in wild orchid populations. It can be found in a wide variety of orchid genera but does not infect plants other than orchids.
--------------	--

**Keywords**

Cymbidium mosaic virus; CymMV; Orchid mosaic virus

---