



Anti-Pspi polyclonal antibody [AP] (CABT-BL6292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Antiserum is developed in rabbit using pure cultures of bacteria as immunogen.
Source/Host	Rabbit
Purification	caprylic acid purified
Conjugate	AP
Applications	ELISA
Format	Liquid
Size	200 µl, 1 ml
Buffer	In 0.05 M TRIS buffer, pH 8.0 and 50% glycerol.
Preservative	None
Storage	The antibody should be stored at 2-8°C. For storage longer than one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-life of 2 years after date of purchase.

BACKGROUND

Introduction	<i>Pseudomonas syringae</i> pv. <i>psis</i> has long been known to be present in pea growing areas and to cause serious crop losses.
Keywords	<i>Pseudomonas syringae</i> pv <i>psis</i> ; Pspi; Pea bacterial blight