



# Anti-Pssy polyclonal antibody [AP] (CABT-BL6294)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Antiserum is developed in rabbit using pure cultures of bacteria as immunogen.
<b>Source/Host</b>	Rabbit
<b>Purification</b>	caprylic acid purified
<b>Conjugate</b>	AP
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	200 µl, 1 ml
<b>Buffer</b>	In 0.05 M TRIS buffer, pH 8.0 and 50% glycerol.
<b>Preservative</b>	None
<b>Storage</b>	The antibody should be stored at 2-8°C. For storage longer then one year, the solution may be frozen at -20°C. Repeated freezing and thawing is not recommended. The solution may be frozen in aliquots if necessary. The antibody has a shelf-live of 2 years after date of purchase.

## BACKGROUND

<b>Introduction</b>	The bacteria <i>Pseudomonas syringae</i> pv. <i>syringae</i> is a member of a large genus of plant-pathogenic bacteria. The bacteria exist as surface dwellers (epiphytes) on leaves and, during wet weather in spring or early summer, can enter through the leaf pores (stomata), causing infections to develop in the young leaves. As the leaf matures these infections cease to expand and are revealed as small patches of dead tissue. As the leaf expands fully, the live tissues pull
---------------------	--

away from the dead patch, which drops out, leaving a 'shothole'.

---

**Keywords**

*Pseudomonas syringae* pv *syringae* (strain specific); Pssy

---