



Anti ABIN1 polyclonal antibody (CABT-BL6311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects ABIN1 at ~72 kDa
Immunogen	human ABIN1 (residues 1-636) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction	This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this
---------------------	---

GENE INFORMATION

Entrez Gene ID [10318](#)

UniProt ID [B7Z8K2](#)
