



Anti-Cullin 5 (aa 577-689) polyclonal antibody (CABT-BL6319)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects Cullin5 at ~91 kDa
Immunogen	human Cullin5 (residues 577-689) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

GENE INFORMATION

Entrez Gene ID [8065](#)

UniProt ID

[Q93034](#)
