



Anti-HOIL1 polyclonal antibody (CABT-BL6324)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects HOIL1 at ~58 kDa
Immunogen	human HOIL1 (residues 1-510) [His-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction	The protein encoded by this gene is similar to mouse UIP28/UbcM4 interacting protein. Alternative splicing has been observed at this locus, resulting in distinct isoforms. [provided by RefSeq, Jul 2008]
---------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

GENE INFORMATION

Entrez Gene ID [10616](#)

UniProt ID [Q9BYM8](#)
