



Anti-HOIP polyclonal antibody (CABT-BL6325)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects HOIP at ~120 kDa
Immunogen	human HOIP (residues 1-1072) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction	The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. [provided by RefSeq, Jul 2008]
--------------	--

GENE INFORMATION

Entrez Gene ID [55072](#)

UniProt ID [Q96EP0](#)
