



## Anti-AMPK alpha 2 (aa 352–366) polyclonal antibody (CABT-BL6337)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	detects AMPK alpha 2 at ~62 kDa
Immunogen	AMPK alpha 2 (residues 352-366) [CMDDSAMHIPPGLKPH]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required