



# Anti-IKK beta polyclonal antibody (CABT-BL6338)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	detects IKK beta at ~87 kDa
<b>Immunogen</b>	human IKK beta (residues 1-736) [6His-tagged]
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	Human
<b>Purification</b>	affinity-purified using immobilized immunogen
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IP
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	phosphate-buffered saline
<b>Preservative</b>	None
<b>Storage</b>	12 months at -20°C; aliquot as required