



Anti-SGK3 polyclonal antibody (CABT-BL6351)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects SGK3 at ~57.1 kDa
Immunogen	human SGK3 (residues 1-300) [GST-tagged]
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction	<p>This gene is a member of the Ser/Thr protein kinase family and encodes a phosphoprotein with a PX (phox homology) domain. The protein phosphorylates several target proteins and has a role in neutral amino acid transport and activation of potassium and chloride channels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.</p>
--------------	--

GENE INFORMATION

Entrez Gene ID [23678](#)

UniProt ID [A0A024R807](#)
