



Anti-UBE2T polyclonal antibody (CABT-BL6355)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	detects Ube2T at ~23kDa
Immunogen	human Ube2T (residues 1-197)
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	affinity-purified using immobilized immunogen
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	None
Storage	12 months at -20°C; aliquot as required

BACKGROUND

Introduction	The covalent conjugation of ubiquitin to proteins regulates diverse cellular pathways and proteins. Ubiquitin is transferred to a target protein through a concerted action of a ubiquitin-
--------------	---

activating enzyme (E1), a ubiquitin-conjugating enzyme (E2), such as UBE2T, and a ubiquitin ligase (E3) (Machida et al., 2006 [PubMed 16916645]).[supplied by OMIM, Mar 2008]

GENE INFORMATION

Entrez Gene ID [29089](#)

UniProt ID [A0A024R9A9](#)
