



# Anti-mono and polyubiquitylated conjugates monoclonal antibody, clone GL3 [HRP] (CABT- BL6369)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Recognises mono and polyubiquitylated conjugates. Does not cross-react with free ubiquitin.
<b>Immunogen</b>	Crude preparation of polyubiquitylated lysozyme
<b>Isotype</b>	IgM
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	GL3
<b>Conjugate</b>	HRP
<b>Applications</b>	WB
<b>Molecular Weight</b>	~150 kDa
<b>Format</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Buffer</b>	10 mM phosphate buffer, 0.15 M NaCl pH 7.4, 0.1% sodium azide
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	12 months at -20°C; aliquot as required