



Anti-Polyubiquitylated conjugates monoclonal antibody. clone GL2 (CABT-BL6370)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognises polyubiquitylated conjugates. Does not cross-react with monoubiquitylated proteins or free ubiquitin.
Immunogen	Crude preparation of polyubiquitylated lysozyme
Isotype	IgM
Source/Host	Mouse
Species Reactivity	Human
Clone	GL2
Conjugate	Unconjugated
Applications	WB
Molecular Weight	~150 kDa
Format	Liquid
Concentration	1mg/ml
Buffer	10 mM phosphate buffer, 0.15 M NaCl pH 7.4, 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	12 months at -20°C; aliquot as required